

Orbicella faveolata and *Siderastrea siderea* cores in nutrients and heat stress

Ana Palacio, Caroline Dennison, Celia Leto, Prati Rosen

Contents

Collection and maintenance	3
Aug 31, 2017	3
Sep 1, 2017	3
Sep 2-12, 2017	4
Sep 18-19, 2017	4
Sep 21, 2017	4
Sep 22, 2017	4
Sep 26, 2017	5
Sep 29, 2017	5
Oct 07, 2017	5
Oct 08, 2017	5
Oct 09, 2017	5
Oct 10, 2017	5
Oct 15, 2017	5
Oct 18, 2017	6
Oct 19, 2017	6
Oct 20, 2017	6
Oct 23, 2017	7
Oct 25, 2017	7
Oct 26, 2017	7
Nov 06, 2017	7
Nov 10, 2017	7
Nov 15, 2017	8
Nov 16, 2017	8
 NUTRIENT ADDITION	 8
Nov 17, 2017	8
Nov 19, 2017	9
Nov 20, 2017	9
Nov 22, 2017	9
Nov 23, 2017	9
Nov 25, 2017	9
Nov 27, 2017	10
Nov 29, 2017	10
Nov 30, 2017	10
Dec 06, 2017	10
Dec 8, 2017	10
Dec 10, 2017	11
Dec 11, 2017	11
Dec 12, 2017	11
Dec 13, 2017	11
Dec 14, 2017	12
Dec 15, 2017	12
Dec 16, 2017	13

Dec 17, 2017	13
Dec 19, 2017	13
Dec 20, 2017	13
Dec 21, 2017	14
Dec 22, 2017	14
Dec 23, 2017	14
Dec 27, 2017	14
Dec 28, 2017	15
Dec 29, 2017	15
Jan 03, 2018	15
Jan 04, 2018	15
Jan 05, 2018	15
Jan 08, 2018	16
Jan 10, 2018	16
Jan 11, 2018	16
Jan 16, 2018	16
Jan 17, 2018	16
Jan 18, 2018	17
Jan 19, 2018	17
Jan 22, 2018	17
Jan 23, 2018	17
Jan 24, 2018	17
Jan 25, 2018	18
Jan 26, 2018	18
Jan 28, 2018	18
Jan 29, 2018	19
Jan 30, 2018	19
Feb 01, 2018	19

RAMP UP 19

Feb 02, 2018	19
Feb 03, 2018	20
Feb 05, 2018	20
Feb 06, 2018	20
Feb 07, 2018	21
Feb 08, 2018	21
Feb 09, 2018	21
Feb 12, 2018	21
Feb 13, 2018	22
Feb 14, 2018	22
Feb 15, 2018	22
Feb 18, 2018	23
Feb 19, 2018	23
Feb 20, 2018	24
Feb 21, 2018	24
Feb 22, 2018	24
Feb 23, 2018	24
Feb 25, 2018	25
Feb 26, 2018	25
Feb 27, 2018	25
Feb 28, 2018	26
Mar 01, 2018	26
Mar 03, 2018	26
Mar 05, 2018	27

Mar 06, 2018	27
Mar 07, 2018	27
Mar 08, 2018	28
Mar 09, 2018	28
Mar 10, 2018	28
Mar 12, 2018	29
Mar 13, 2018	29
Mar 14, 2018	29
Mar 15, 2018	30
Mar 19, 2018	30
Mar 20, 2018	30
Mar 21, 2018	30
Mar 22, 2018	31
Mar 23, 2018	31
Mar 26, 2018	31
Mar 27, 2018	32
Mar 28, 2018	32
April 4, 2018	32
April 5, 2018	32
Apr 12, 2018	33
Apr 13, 2018 (Friday the 13th)	33
Apr 20, 2018	33
Apr 26, 2018	33
Apr 30, 2018	34
May 07, 2018	34
May 16, 2018	34
May 17, 2018	34
May 21, 2018	35
May 22, 2018	35

Collection and maintenance

Aug 31, 2017

Cleaned and organized all the Ssid cores

- Master file updated with the cores in the tank
-

Sep 1, 2017

Ssid IPAM T= 0 and Pictures with color scale

- Colony 20 by Becca and AP
- Picture # 2 has more cores than the spread sheet, check numbers :/
- Colony 22, 23, 24, 27, 28 and 30 by PR and AP

Ssid Tissue samples T= 0

- Selected cores from colony 20 by Becca and AP

- Selected cores from colony 22, 23, 24, 28 and by PR and AP
- Colonies 27 and 30 were NOT sampled!!

Ssid SDS incubation

- All the samples were incubated in 1% SDS for 1:30 hours
-

Sep 2-12, 2017

IRMA!!

- The corals were not cleaned or fed until 13th when AB and Maddie could get into RSMAS and clean them
-

Sep 18-19, 2017

- Cleaned and fed
 - Lots of mortality
-

Sep 21, 2017

Ssid Tissue samples T= 0 for colonies Ssid 27

- Tissue samples taken by Carly Dennison
 - SDS incubation
-

Sep 22, 2017

Ssid Tissue samples T= 0 for colonies Ssid 30

- Tissue samples taken by CD
- SDS incubation

Ofav new cores from colony Ofav “33”

- O.fav 33
 - This colony has been kept in the wet lab for a long time, one side is pale/bleach since It never got light
 - The side facing up had some tissue loss beacause it was out of the water
-

Sep 26, 2017

Ofav 33 colony moved to tank 5A

- Colony cleaned with Bayer and moved to tank 5A
-

Sep 29, 2017

Ssid cores (8 per genotype) were selected to be used in the nutrient experiment

- Selection of cores
- Made new 7*8 eggcrates for them ———

Oct 07, 2017

- Make stock solutions and started the pums.
 - Took water samples for nutrient analysis day 1
-

Oct 08, 2017

- Took water samples for nutrient analysis day 2
-

Oct 09, 2017

- Took water samples for nutrient analysis day 3
 - Cleaned all Ssid cores
-

Oct 10, 2017

IPAM 1 for Ssid and Pictures w color scale

- All the Ssid cores were ok
 - Values are lower vs Sep 1st (IPAM 0)
-

Oct 15, 2017

New Nutrient stocks

- Made more Stock 1 and Stock 2 Nutrients
 - Stock 1 N: 4.011825g of NH_4Cl + 250mL of RO water
 - Stock 1 P: 1.034925g of $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ + 250mL of RO water

- Stock 2 N: 21mL of Stock 1 + 987mL of RO water
 - Stock 2 NP: 21mL of N_Stock 1 + 42mL of P_Stock 1 + 945mL of RO water
 - New Stock 2 solutions were attached to the pumps at 4:06PM to follow volume delivery for 24hours
-

Oct 18, 2017

IPAM 2a for Ssid

- Data taken by AP
- Values are lower in the replicate 2 vs 1! , do not look good to start the experiment

Tissue samples 2 for Ssid

- Samples taken by AP
- SDS incubation
- Core 28-19 was not sampled... looks dead?

Buoyant weight 2 for Ssid

- Replicate 1 by CD and PR
 - Replicate 2 by CD and AP
-

Oct 19, 2017

IPAM 2b for Ssid

- Data was taken again since yesterdas values were so low
- Less standard variation this time, but mean values did not improve

ProK incubation samples Time 2

- Prok by CD
-

Oct 20, 2017

Organic Extraction

- All Ssid R1 and almost all R2
 - Started by BG, finished by PR
-

Oct 23, 2017

2nd EtOH

- All Ssid R1 and almost all R2

Corals moved inside aquaria, not nutrients added yet

Oct 25, 2017

EtOH wash -> DNA/0.1X TE

- Ssid R1 (except the first set of 48) and almost all R2, by BG and CD

Organic Extraction last Ssid R2 and Acer R1

- CTAB CD, Chl PR, 1st ETOH CD
 - Last 19 samples of Ssid R2
-

Oct 26, 2017

EtOH wash -> DNA/0.1X TE

- First Ssid R1 (48 samples) by AP

2nd EtOH last 19 samples in Ssid R2

- Used the NaOH made yesterday and added a - control to the extraction
-

Nov 06, 2017

IPAM 3 for Ssid

- Values look better!
 - min=0.47 (22-24), mean and median ~0.55, max 0.59 (28-23)
-

Nov 10, 2017

Buoyant weight before sampling for Ssid R1 (Control and N)

- Data taken by CD and BG
-

Nov 15, 2017

Buoyant weight before sampling for Ssid R1 (NP) and R2

- R1 NP by BG and PR
- R2 Control by BG and PR
- R2 N by PR
- R2 NP by CD

Tissue samples T=4

- Ssid samples taken by CD
- Ofav samples taken by PR

Pictures T=4 for Ssid

Nov 16, 2017

IPAM T=4 for Ssid and Ofav

- Ssid values declined, the increased in T=3 is reverted and values look similar to T=2
 - It was maybe sampling/handeling?
-

NUTRIENT ADDITION

Nov 17, 2017

STARTED NUTRIENT ADDITION!!!

- To start concentrations:
 - N: 1mL of NH_4Cl [Stock 1 = 300mM] in 30L filtered sea water
 - NP: 1mL of NH_4Cl [Stock 1 = 300mM] + 1mL NaH_2PO_4 [Stock 1 = 30mM] in 30L filtered sea water
- To keep concentrations
 - N: 24mL of NH_4Cl / day [Stock 2 = 6.25uM]
 - NP: 24mL of NH_4Cl + NaH_2PO_4 / day [Stock 2 = 6.25uM + 1.25uM]

Split Ofav cores

- Cores were split to balance YII values in each treatment
 - Some were not labeled, so randomly allocated
-

Nov 19, 2017

Fed Ssid and ofav

- Outside the aquariums with the larvae food + reef chilli

Insecticide dip for Ofav and Revive dip for Ssid

Nov 20, 2017

Ssid ProK incubation

- All Ssid samples by CD

Ofav moved inside the aquariums

- Check labels and re-label some
 - Pictures with color scale
-

Nov 22, 2017

Fed Ssid and Ofav

- Inside the aquariums with reef chilli
-

Nov 23, 2017

IPAM time = 5

- AP and CD
- Values in Ofav colony 34 declined
- Values in Ssid increased, more in R2

Water change and horizontal rotation

Nov 25, 2017

Initial buoyant weight Ofav T= 5

- After initial samples and already in Nutrients
-

Nov 27, 2017

ProK incubation o all Ofav samples (initial time point)

- 100uL aliquots + 6uL of ProK by RG and AP
-

Nov 29, 2017

IPAM time = 6

- Values in Ofav increased (AP and BG)
- Values in Ssid increased back to time 4! (AP and PR)

Water change and horizontal rotation

Water sample for nutrient analysis

* ~ 50 mL from each aquarium

Nov 30, 2017

Measured light in all aquaria

Water samples day 2 for nutrient analysis

* ~ 50 mL from each aquarium

Dec 06, 2017

IPAM time = 7

- PR and CD
- ~ Constant values

Water change and vertical rotation

Dec 8, 2017

CTAB first batch of Ssid sampled in 11-15-2017

- 1st block Ssid 20-1 to 28-3
- 2nd block Ssid 28-6 to 24-23

- No negative control, started by BG
-

Dec 10, 2017

1st ETOH first batch of Ssid sampled in 11-15-2017

- AP
-

Dec 11, 2017

2nd ETOH first batch of Ssid sampled in 11-15-2017

- AP
-

Dec 12, 2017

ETOH wash first batch of Ssid sampled in 11-15-2017

- AP
-

Dec 13, 2017

1st ETOH

- second batch of Ssid (65 samples) + first batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-1 to 34-30s)
- CTAB by CD
- We ran these O.fav samples and they have late coral CTs, we do not know why!!!!

IPAM time = 8

- CD and AP
- Ssid values increased a bit, specially for N, Control looks lower
- Ofav values increased a bit, Control looks lower

Water change and horizontal rotation

Dec 14, 2017

2nd ETOH:

second batch of Ssid (65 samples) + first batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-1 to 34-30)
by CD

Organic Extraction and 1st ETOH

- second batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-1 to 34-49, Nutrients and CO2) + Ofav CO2
- by CD

Ofav BW before samples

- Data taken by CD and my mom

Ssid tissue samples + Picture with color scale

- Taken by CD
-

Dec 15, 2017

ETOH wash

- second batch of Ssid (65 samples) + first batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-1 to 34-30)
- by AP

2ndt ETOH

- second batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-30 to 34-49, Nutrients and CO2) + Ofav CO2
- by AP

Ofav tissue samples (T=8) + Picture with color scale

- Taken by PR

Sid BW after samples

- Data taken by CD and my mom

ProK for all Ofav and all Ssid R1 up to 30-14

Dec 16, 2017

ETOH wash

- second batch of Ofav sampled in 11-15-2017 (Colony 34 from 34-30 to 34-49, Nutrients and CO2) + mess of Ofav in CO2 and Nutrients
- by AP

Feeding and cleaning

- All spp fed with green food.
-

Dec 17, 2017

Water change and vertical rotation

- AP

New Stock 2 of nutrients

- AP
-

Dec 19, 2017

ProK incubation Ssid R2 from Dec 14

- from 20-26 to 20-30 (55 samples) with last Acer samples
-

Dec 20, 2017

Organic extraction Ssid R2 from Dec 14

- from 20-26 to 20-30 (55 samples) with last Acer samples
- CTAB AP- 1st EtOH PR

ProK incubation last of Ssid R2 from Dec 14

- from 20-36 to 30-28 (26 samples)
 - PR
-

Dec 21, 2017

Second EtOH of Ssid R2 from Dec 14

- from 20-26 to 20-30 (55 samples) with last Acer samples
- AP

Cleaned tanks

- AP, corals left without nutrients
-

Dec 22, 2017

EtOH wash of Ssid R2 from Dec 14

- from 20-26 to 20-30 (55 samples) with last Acer samples
- AP

Cleaned corals + horizontal rotation

- AP
-

Dec 23, 2017

Organic extraction of all Ofav samples from Dec 14 and some R2 Ssid

- All Ofav samples
 - Ssid from 20-36 to 28-23 (21 samples)
-

Dec 27, 2017

Second EtOH Ofav samples from Dec 14 and some R2 Ssid

- All Ofav samples
- Ssid from 20-36 to 28-23 (21 samples)
- Worked with S.Diaz de Villegas

Fed corals + vertical rotation

Dec 28, 2017

Organic Extraction and 1st EtOH R1 Ssid and 4 last samples R1 (from Dec 14)

- AP and NC
-

Dec 29, 2017

2nd EtOH R1 Ssid and 4 last samples R1 (from Dec 14)

- AP and NC
-

Jan 03, 2018

Fed corals + horizontal rotation

- With S.D

IPAM time = 9

- SD and AP
 - All values declined for Ofav, except N R2
 - All values declined for Ssid.
-

Jan 04, 2018

EtOH wash R1 Ssid and 4 last samples R1 (from Dec 14)

- Worked with S.Diaz de Villegas
 - DNA in 0.1 TE
-

Jan 05, 2018

New nutrient stock (1) and (2)

- Replaced stock 2 in pumps
-

Jan 08, 2018

Water change with non filtered water

- Added more stock 1
-

Jan 10, 2018

Fed corals and stayed in extra tank for cleaning next day

- CD and SD
-

Jan 11, 2018

Cleaned corals + vertical rotation

- CD and SD
-

Jan 16, 2018

Cleaned tanks + horizontal rotation

- CD and SD
-

Jan 17, 2018

qPCR to check A and B presence

- Plate 1 Ofav (StepOnePlus)
 - Plate 1 (New machine) and 2 Ssid (StepOnePlus)
 - Ran by CD with SYBR (expired in 2010) instead TaqMan
 - All the plates show amplification of A and B, including in the NTC, I think the primers create primer dimer and the SYBR was not specific enough showing increasing florescence signal.
 - Melting curved showed multiple and and lower temperature peaks than the positive controls.
 - Some of the Ofav samples showed earlier curves for B, re-run these plate with TaqMan MM instead SYBR
-

Jan 18, 2018

qPCR to check for A and B using TaqMan MM

- Rerun of Plate 1 (Ofav) using TaqMan MM instead of SYBR
- Same setup
- Run by CD
- Samples from colony 31 and colony 6 have B symbionts. Other colonies are D dominated

Water change

- CD, the corals were left in holding tank to IPAM next day
 - Corals were fed
-

Jan 19, 2018

IPAM T=10 + Pictures

- CD and AP
- Ssid declined, but similar to IPAM before the last one

Corals back to treatments

- Corals moved inside the tanks and nutrients added
 - Vertical rotation
-

Jan 22, 2018

Buoyant weight for Ofav

- CD
-

Jan 23, 2018

Water change by CD

Jan 24, 2018

qPCR Ssid from Dec 14

- Plates 3, 4 and 5 ran by CD in StepOnePlus
- Ssid ran with SYBR and CD with TaqMan MM

- Old + control for plates3, made a new one from the rest of the plates.
-

Jan 25, 2018

qPCR Ssid from Dec 14

- Plates 6 and 7 ran by AP in StepOnePlus

IPAM T=11

- Ofav R2 declined, R1 stable

Pics T=11

- SD
-

Jan 26, 2018

qPCR Ssid from Dec 14

- Plates 7 and 8 ran by AP in StepOnePlus

Nutrient Stock Replenished/HOBOs checked

- Changed the bottles, but used stock made before.

Ofav Tissue samples (Time=12)

- PR
-

Jan 28, 2018

Water change + vertical rotation

- Fed corals
 - Corals (both replicates) had mesenterial filaments out
 - AP and CL
-

Jan 29, 2018

ProK incubation Ofav samples from Jan 26

- By SD and CD
-

Jan 30, 2018

IPAM T=12

- CD and AP

qPCR Ofav samples with no B

- Plates 2 and 3, coral with SYBER, melt curve added.
 - CD
-

Feb 01, 2018

Buoyant weight before samples for Ssid

- CD

Tissue samples for Ssid

- CD
-

RAMP UP

Feb 02, 2018

ProK incubation for Feb 1 Ssid samples

- CD

Ctab and organic extraction of all Ofav samples from Jan 26

- CD and PR

Thermal stress begins

- Temperatures increased 1 degree (26 → 27)
 - AP
-

Feb 03, 2018

Organic and 1st EtOH precipitation for Ssid samples taken on Feb 1

- SD
- 30-9, 30-13, 20-38, 20-45, and 23-54 were all messed up need to be reextracted

2nd EtOH precipitation, EtOH wash, and resuspension for Ofav samples taken on Jan 26

- RESUSPENDED IN 1xTE!!!!
 - SD
-

Feb 05, 2018

Took samples for blasting and put them in -80

- chose 1 sample per genotype (where available)

Water changes and horizontal AND vertical rotation

- done by AP and PR

Took water samples after nutrients were added for analysis

- done by PR

Temperature increase from 27 to 28

- corals were kept at 27 over the weekend

2nd EtOH precipitation for Ssid samples minus the 6 that were messed up over the weekend

- SD
-

Feb 06, 2018

Temperature increase 29

Feb 07, 2018

IPAM time=13

- IPAM done by CD and AP

Water samples taken

- CD

Temperature increase 30

- Corals were at 29 degrees and temperature increased to 30 at the end of the day
-

Feb 08, 2018

Water samples taken

- SD
-

Feb 09, 2018

Buoyant weight after samples Ssid

- NP R2; C, N and NP R2
 - CD and SD
-

Feb 12, 2018

Begin tissue blasting for Ssid samples sacrificed on Feb 5

- 10 samples taken from the -80
- Took aliquots for:
 - Cell counts (1mL + 50uL lugols)
 - DNA (1mL + 10% SDS)
 - Chlorophyll A (~ 2-3mL filtered + 2mL of Methanol)
 - Lipids (~ 2-3mL filtered + moved back to -80)

Salinity and Temperature

- Control 1: 40 - 30.6
- Control 2: 39.9 - 30.6
- N1: 40.3 - 30.7
- N2: 39.8 - 30.6
- NP1: 39.9 - 30.7

- NP2: 35.5 - 30.5

IPAM and Pictures time=14

- IPAM done by CD and AP

Water change

- Vertical and horizontal rotation

Feeding

- PR
-

Feb 13, 2018

Salinity and Temperature

- Control 1: 34.4 - 30.6
 - Control 2: 34.5 - 30.7
 - N1: 34.6 - 30.6
 - N2: 34.4 - 30.7
 - NP1: 34.5 - 30.6
 - NP2: 34.6 - 30.8
-

Feb 14, 2018

Salinity and Temperature

- Control 1: 34.8 - 31.4
- Control 2: 35.4 - 31.7
- N1: 35.3 - 31.4
- N2: 35.2 - 31.6
- NP1: 35.2 - 31.5
- NP2: 35.3 - 31.7

Pictures with colorscale taken of Ofav

Feb 15, 2018

Salinity and Temperature

- Control 1: 36.6 - 30.6

- Control 2: 36.7 - 30.6
- N1: 36.4 - 30.9
- N2: 36.2 - 30.6
- NP1: 36.3 - 30.9
- NP2: 36.5 - 30.7

IPAM and pictures with color scale time=16

- IPAM done by CD and AP

Water changes

- Vertical and horizontal shift of aquaria
- Done by AP and CD

Stopped nutrient addition

Feb 18, 2018

Added fresh water to aquaria

~ 100mL of RO water was added to each one

Feb 19, 2018

Salinity and Temperature

- Control 1: 37.4 - 31.1
- Control 2: 38.1 - 31.5
- N1: 37.9 - 31.2
- N2: 38.2 - 31.5
- NP1: 38.2 - 31.2
- NP2: 38.0 - 31.5

IPAM and pictures with color scale time=17

- IPAM done by CD and AP

Water changes + Feeding

- Vertical and horizontal shift of aquaria
 - Done by AP and PR
-

Feb 20, 2018

Salinity and Temperature

- Control 1: 35.1 - 31.5
 - Control 2: 35.5 - 31.8
 - N1: 35.3 - 31.6
 - N2: 35.4 - 31.8
 - NP1: 35.5 - 31.6
 - NP2: 35.4 - 31.7
-

Feb 21, 2018

Salinity and Temperature

- Control 1: 36.1 - 31.5
 - Control 2: 36.0 - 31.7
 - N1: 36.0 - 31.7
 - N2: 36.1 - 31.8
 - NP1: 36.2 - 31.6
 - NP2: 36.1 - 31.7
 - RO water added to aquaria to reach target salinity (34.5-35.0)
-

Feb 22, 2018

IPAM and pictures with color scale time=18

- IPAM done by CD and AP

Water changes and feeding

Got HOBOS data from Jan-Feb

Feb 23, 2018

Razor blade samples taken

- Ssid done by CD
- Ofave done by Prati

Salinity and Temperature

- Control 1: 35.4 - 31.7
 - Control 2: 35.7 - 31.6
 - N1: 35.6 - 31.8
 - N2: 35.8 - 31.6
 - NP1: 35.9 - 31.8
 - NP2: 35.8 - 31.6
 - RO water added to aquaria to reach target salinity (34.0-34.3)
-

Feb 25, 2018

- Added ~ 1 cup of RO water
-

Feb 26, 2018

Salinity and Temperature

- Control 1: 35.8 - 31.7
- Control 2: 35.8 - 31.6
- N1: 36.0 - 31.7
- N2: 36.2 - 31.6
- NP1: 36.4 - 31.7
- NP2: 36.3 - 31.6

IPAM and pictures with color scale time=18

- done by CD, AP, and SD
-

Feb 27, 2018

Salinity and Temperature

- Control 1: 34.4 - 31.5
- Control 2: 34.6 - 31.6
- N1: 34.9 - 31.6
- N2: 34.8 - 31.7
- NP1: 35.0 - 31.6
- NP2: 34.8 - 31.7

CTAB incubation, organic extraction, and 1st EtOH precipitation

- Ofav and Ssid samples taken on Feb 23, 2018 (Done by CD and SD)
-

Feb 28, 2018

Salinity and Temperature

- Control 1 35.5 - 31.6
- Control 2: 35.6 - 31.6
- N1: 35.8 - 31.6
- N2: 35.7 - 31.6
- NP1: 35.8 - 31.6
- NP2: 35.8 - 31.6

2nd EtOH precipitation for all Ofav and Ssid samples

- Done by AP and CD
-

Mar 01, 2018

Water changes and feeding

- Vertical and horizontal rotation

IPAM and Pics T=19

- Done by CD and AP

EtOH wash for samples taken on Feb 23, 2018

- Done by CD and SD
-

Mar 03, 2018

Salinity and Temperature (adjusted with ~ 2 cups of RO)

- Control 1 36.2 - 31.3 (35.3)
- Control 2: 34.4 - 31.6 (35.4)
- N1: 35.8 - 31.7 (35.4)
- N2: 36.2 - 31.6 (35.4)
- NP1: 36.4 - 31.7 (35.3)
- NP2: 36.3 - 31.6 (35.4)

ProK incubation dying corals

Mar 05, 2018

Blasted The rest of Ssid fragments

- done by CD, AP, SD and CL

Water changes and feeding

- Vertical and horizontal rotation

IPAM and Pics T=20

- Done by CD and AP

SDS and Prok incubation of blastate-samples

Mar 06, 2018

Chlorophyll for all Ssid fragments blasted on Mar 05, 2018

Tissue samples for Ofav and Ssid

- CD took Ssid and SD took Ofav

qPCR Ofav with BCD

- Plates # ??? samples from date???
-

Mar 07, 2018

Salinity and Temperature

- Control 1: 37.0 - 31.5
- Control 2: 37.0 - 31.7
- N1: 37.0 - 31.5
- N2: 37.2 - 31.7
- NP1: 37.1 - 31.6
- NP2: 37.1 - 31.8

qPCR run for Ofav colonies 6 and 31 from December sample date

- Run by CD

ProK incubation 3/6/18 samples

Mar 08, 2018

Female corals strike

Salinity and Temperature

- Control 1: 35.9 - 31.4
- Control 2: 35.4 - 31.6
- N1: 35.6 - 31.3
- N2: 35.7 - 31.6
- NP1: 35.5 - 31.4
- NP2: 35.7 - 31.7

IPAM T=20?

- AP and CD

Took fragments to blast and put them in -80

- AP and CD
-

Mar 09, 2018

qPCR for Ssid samples

- 27-87 to 30-38 from 12/14/17 the rest of the samples are from 2/1/17

CTAB incubation, organic extraction, and 1st EtOH precipitation

- Ssid samples up to 30-21 (done by CD and PR)
-

Mar 10, 2018

4 qPCR plates for Ssid samples

- 2/1/18
- Ssid Pax-C sybr green, C/D TaqMan

- Step One machine
 - Run by CD
-

Mar 12, 2018

2 qPcr plates for Ssid samples

- 2/1/18
- Ssid Pax-C sybr green, C/D TaqMan
- Step One machine
- Run by CD

Rerun of Ofav Plate 5 in qPCR

- original plate had amplification in NTC ** Primer dimer
 - Step One machine
 - 0.02 threshold with sybr
 - Run by CD
-

Mar 13, 2018

CTAB incubation, Organic extraction, and 1st EtOH for Ofav and second part of Ssid samples taken on 3/6/18

- Done by CD
-

Mar 14, 2018

Tissue blasting for all Ofav samples sacraficed on Feb 5

- Samples taken from the -80
- Took aliquots for:
 - Cell counts (1mL + 50uL lugols)
 - DNA (1mL + 500DNAB + 83uL of 10% SDS) SDS incubation for 1:30 at 65C
 - Chlorophyll A (~ 2-3mL filtered + 4mL of Methanol)
 - Lipids (~ 2-4mL filtered + moved back to -80)

2nd ETOH for Ofav and last Ssid samples taken on 3/6/18

- Done by SD and CD
-

Mar 15, 2018

Chl-a for all Ofav samples sacrificed on Feb 5 (before bleaching)

Corals in recovery moved outside the aquaria and put on recirculating system

- Ssid and Ofav dip in insecticide

Lights with calibrated lightmeter

- Control 1: Right top: 321
 - Control 2: Right bottom: 420
 - N1: Left, top: 350
 - N2: Left, bottom: 365
 - NP1: Center, top: 345
 - NP2: Center, bottom: 435
-

Mar 19, 2018

EtOH wash and resuspension of DNA samples taken on 3/6/18

- Done by CD
-

Mar 20, 2018

Blasted all Ofav fragments collected on Mar 08 (bleaching)

- done by CD, AP, SD
- Samples taken from the -80
- Took aliquots for:
 - Cell counts (0.5mL + 50uL lugols)
 - DNA (0.5mL + 1/2(500DNAB + 83uL of 10% SDS)) SDS incubation for 1:30 at 65C
 - Chlorophyll A (~ 3-4mL filtered + 4mL of Methanol)
 - Lipids (~ 3-4mL filtered + moved back to -80), very few samples had enough blastate for lipids

SDS incubation of blastate samples

Mar 21, 2018

Chl-a for Ofav fragments collected on Mar 08

- done by CD, AP

PCR ITS-2 selected samples to sequence symbionts

- fragments per genotype
- Done by SD

Cell counts for before bleaching

- Ssid by CL
-

Mar 22, 2018

qPCR for Ssid in Nutrients Plate 17 and 18

Re-extract Ofav 20-29 (Pro-K) because no DNA sample was found

- ProK with blastate samples
- Done by SD

ProK incubation of Ofav blastate taken on Mar 08, 2018 (bleaching)

- Done by SD

PCR ITS-2 gels to sequence symbionts

- Done by SD
-

Mar 23, 2018

qPCR for Ssid in Nutrients Plates 19 and 20

- Done by CD
-

Mar 26, 2018

Blasted all the Ssid cores after bleaching

- 42 cores
- CL, CD and AP
- Samples taken from the -80
- Took aliquots for:
 - Cell counts (0.5mL + 50uL lugols)
 - DNA (0.5mL + 1/2(500DNAB + 83uL of 10% SDS)) SDS incubation for 1:30 at 65C
 - Chlorophyll A (~ 3-6mL filtered + 4mL of Methanol)
 - Lipids: only 2 samples had enough blastate for lipids

SDS incubation of blastate samples

Mar 27, 2018

Chl-a extraction for Ssid cores after bleaching

- 42 cores
- SD, CD and AP

ProK incubation Ssid blasted on Mar 26

- SD

IPAM Ofav and Ssid after 3 weeks of recovery

- AP and CD
-

Mar 28, 2018

Organic extraction for Ssid and Ofav blasted cores after bleaching

- SD
-

April 4, 2018

qPCR plate 22 for Ssid Nutrients

- Checking if amplification cycles are consistent from time point to time point
 - Chose 1 sample per genotype and ran it from each time point to compare
 - CD
-

April 5, 2018

qPCR plates 6-9 for Ofav Nutrients

- Ofav sybr green SC assay
 - Taqman assay for C/D and B (plates 8 and 9)
-

Apr 12, 2018

Water samples taken from experiment tank and cleaning tank (5B) to measure DIC and TA for insight to what CO₂ levels are in system

- Samples taken from incoming water, recirculating system, aquarium with pumps, and aquarium with pump and bubbler
- 8 samples taken total (1 replication)
- Done by Emma Pontes and new lab manager in Langdon Lab

2nd EtOH precipitation done for one set of Ssid

- CD

qPCR for Ssid Plate 24

- Vac fug broken, hand spin *CD
-

Apr 13, 2018 (Friday the 13th)

qPCR for Ssid Nutrients plates 25, 27-31

- Ssid syrb green assay
- Taqman for C/D
- Done by CD

2nd EtOH precipitation for the remaining Ofav and Ssid sets

- CD
-

Apr 20, 2018

IPAM (T23) Ofav and Ssid after ~6 weeks of recovery

- AP and CD
-

Apr 26, 2018

1:2 dilutions for 4 cell count samples to recount

- 34-38, 6-86, 20-21, and 6-118
- high SE and >50 cells per chamber
- 20uL of sample, 20uL DNAB
- done by CD

Apr 30, 2018

qPCR for Ssid Nutrients Plates 32-35

- Done by CD
-

May 07, 2018

qPCR for Ofav Nutrients Plates 10-16

- Done by CD
-

May 16, 2018

New positive (+) control made for qPCR with A,B,C,D symbionts

- 20 uL of OfavSC (Ofav Plate 16 columns 7-12) qPCR product
 - 10 uL of A1 standard actin cleaned PCR product
 - 10 uL of A2 standard actin cleaned PCR product
 - 10 uL of B actin standard #3 PCR product
 - 10 uL of B actin standard #4 PCR product
 - 10 uL of D actin standard #3 PCR product
 - 10 uL of D actin standard #4 PCR product
 - 20 uL of Acer CaM (Acer Plate 31 rows BCD column 1)
 - 10 uL of C standard #3 M13 PCR product
 - 10 uL of C standard #4 M13 PCR product
 - Total volume of 100 uL then cleaned using Wizard SV Gel and Clean-up System protocols
 - Done by CD and AP
-

May 17, 2018

qPCR for Ofav Nutrients from 03/06/18

- Check for A,B,C,D symbionts post bleaching and to rule out contamination
 - Master mix: Actin primers and probes and EMM
-

May 21, 2018

PCR stock made

- Done by CD

1:10 dilutions of the super positive made on 05/16/18

- 10 uL of (+) and 90 uL of UP water
-

May 22, 2018

PCR for Ssid, Ofav, and Acer Nutrient using Ssid PaxC, OfacSC and AcerCaM primers respectively

- this will be used for copy #, fluormetry corrections, and specialized/species super positive qPCR sample
- Done by CD

qPCR for Ssid Nutrients from 03/06/18

- Check for A,B,C,D symbionts post bleachin and to rule out contamination
 - Master mix: Actin primers an probes and EMM
 - 1:10 super (+) used
 - Done by CD
-